

DESIGNED FOR SAFETY. DESIGNED TO FIT.


Shiley™ XLT extended-length
tracheostomy tubes



Proximal
extension



Distal
extension



Medtronic
Further. Together

SIZING AND ORDERING INFORMATION

For Shiley™ XLT extended-length tracheostomy tubes

DESCRIPTION							CATALOG NUMBER		
Extension	I.D.	O.D.	Total length	Proximal length	Radial length	Distal length	Cuffed	Cuffless	Disposable inner cannula
Distal	5.0 mm	9.6 mm	90 mm	5.0 mm	37.0 mm	48.0 mm	50XLTCD	50XLTUD	50XLTIN
Proximal	5.0 mm	9.6 mm	90 mm	20.0 mm	37.0 mm	33.0 mm	50XLTCP	50XLTUP	50XLTIN
Distal	6.0 mm	11.0 mm	95 mm	8.0 mm	38.0 mm	49.0 mm	60XLTCD	60XLTUD	60XLTIN
Proximal	6.0 mm	11.0 mm	95 mm	23.0 mm	38.0 mm	34.0 mm	60XLTCP	60XLTUP	60XLTIN
Distal	7.0 mm	12.3 mm	100 mm	12.0 mm	39.0 mm	49.0 mm	70XLTCD	70XLTUD	70XLTIN
Proximal	7.0 mm	12.3 mm	100 mm	27.0 mm	39.0 mm	34.0 mm	70XLTCP	70XLTUP	70XLTIN
Distal	8.0 mm	13.3 mm	105 mm	15.0 mm	40.0 mm	50.0 mm	80XLTCD	80XLTUD	80XLTIN
Proximal	8.0 mm	13.3 mm	105 mm	30.0 mm	40.0 mm	35.0 mm	80XLTCP	80XLTUP	80XLTIN

Tracheostomy tubes 1/box, disposable inner cannula 10/box. Supplied with both cotton twill tape and adjustable tracheostomy tube holder (TTH).

© 2017 Medtronic. All rights reserved. Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company. 01/2017-16-AW-0079-[WF#1419324]

6135 Gunbarrel Avenue
Boulder, CO 80301

800.635.5267

medtronic.com/covidien

Medtronic